



**-25°C FREEZER CHEST BFRZ-25-404**

## -25°C FREEZER CHEST BFRZ-25-404

A perfect partner for your multi cold storage purposes in a contemporary lab set up. Microprocessor based temperature controls provide outstanding temperature uniformity with high quality compressor for energy efficiency. Safety features and durable interiors provide you with the most dependable performance.

Used in Laboratories, Routine Cold Storage, Biologicals, Pharmaceuticals, Clinical, Research.

Also known as Laboratory ULT freezer, ULT Freezer, Ultra-low Freezer.

## BFRZ-25-404 -25°C FREEZER CHEST



Double Door Drawer x 12 USB port

## SPECIFICATIONS

Model	BFRZ-25-404
Capacity	358 L
Temperature Range	-10°C to -25°C
Alarm	High&low temp, power failure, control error, door open, high ambient temperature
External Dimension	720×830×1870 mm
Power	187~242V, 50/60Hz
Optional Parts	110V transformer; Datalogger



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)