



STANDARD FUME HOOD BFSD-207

STANDARD FUME HOOD BFSD-207

Designed to provide personnel protection from toxic fumes and hazardous volatile vapours, our standard fume hoods are reliable and easy to use. Its inbuilt ventilation system with a front window eliminates fire and explosion risk. Acid\alkali resistant interiors and UV lamp to sterilise the work area, this product is the best choice for reducing potential exposure hazard to the user.

Used in Metal analysis, Pharmacy, Research Field, Commercial, Medical, Laboratory.

Also known as Standard Fume Hood, Laboratory Fume Hood, Chemical Fume Hood.

BFSD-207 STANDARD FUME HOOD



Interior body includes melamine board with acid and alkali resistance function

LCD display

Exterior body is of cold-roll steel with anti-bacterial powder

Working chamber is made of solid chemical-resistant material

UV lamp to sterilize the working area

Adjustable air velocity with electrical controlled glass

Adjustable blower speed of a centrifugal fan

SPECIFICATIONS

Model	BFSD-207
Air Velocity	0.3~0.8 m/s
Max Opening	520 mm
Work Surface Height	750 mm
System Exhaust Volume	990m ³ /h
Front Window	Motorized, 5 mm toughened glass; Anti-UV Exterior Body: cold rolled steel with anti-bacteria powder coating.
Work Table Material	Chemical resistant Phenolic resin
LED Lamp	16Wx 1
UV Lamp	20 W Emission of 253.7 nanometers for most efficient decontamination
Blower	Built-in centrifugal blower, Speed adjustable with 9 levels
Standard Accessory	Water tap, Gas tap, Water sink, base stand, UV lampx2, LED lamp, Socket, 4.3 meters exhaust duct, Active carbon filter, foot switch, Waterproof Socket- 2 pcs
Optional Accessory	HEPA filter
Internal Size (WxDxH)	1380x730x745 mm
External Size (WxDxH)	1500x840x2150 mm
Package (Main Body)	1650x1120x1680 mm
Noise Level	≤60 dB
Gross Weight	290 kg
Consumption	500 W
Power Supply	AC220V±10%, 50/60Hz 110V±10%, 60Hz

ACCESSORIES

Accessory Code	Name
6000610006	Fluorescent lamp
6000610007	UV lamp
6000610008	Base stand
6000610009	Waterproof socket
6000606006	water tap

Accessory Code	Name
6000606007	Gas tap
6000606008	Water sink
6000610010	4 meters PVC exhaust duct
6000606012	Active carbon filter
6000610011	Foot switch



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com