



## UPRIGHT FREEZER DOUBLE DOOR BFUD-25-201

## UPRIGHT FREEZER DOUBLE DOOR BFUD-25-201

High performance refrigerators featuring microprocessor control for temperature accuracy and ergonomic design for easy usage. Designed for space saving storage, it provides outstanding thermal performance along with safety features.

Used in Protection, Frozen Blood, Vaccines, Enzymes, Reagents, Research.

Also known as Double door refrigerated Freezer, Laboratory Combined Refrigerator and Freezer , Laboratory double door refrigerated freezer.

## BFUD-25-201 UPRIGHT FREEZER DOUBLE DOOR

- Microprocessor temperature controller
- LCD display
- Vertical cabinet with two chamber rooms
- ABS liner interior and painted steel board exterior
- 12 drawers made of ABS material inside
- Door with safety lock and rotatable sure-grip door handle
- Polycarbonate insulated inner doors
- Insulation thickness: 60 mm
- Password protected configuration page
- Refrigerant: CFC free R134a
- Danfoss compressor (2 units) and EBM fan motor
- 4 casters for easy handling
- Safety functions: High or low temperature alarm, system failure alarm

## SPECIFICATIONS

Model	BFUD-25-201
Capacity	450 L
Temperature Range	-10°C to -25°C
Overall Dimension	660x780x1850 mm
Package Size	765x850x2050 mm
Weight (Net/Gross)	120/140 kg
Average Power	450 W
Power Supply	220V,50/ 60Hz, 110V,50/60Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)