



UPRIGHT FREEZER DOUBLE DOOR BFUD-40-101

UPRIGHT FREEZER DOUBLE DOOR BFUD-40-101

High performance refrigerators featuring microprocessor control for temperature accuracy and ergonomic design for easy usage. Designed for space saving storage, it provides outstanding thermal performance along with safety features.

Used in Protection, Frozen Blood, Vaccines, Enzymes, Reagents, Research.

Also known as Double door refrigerated Freezer, Laboratory Combined Refrigerator and Freezer , Laboratory double door refrigerated freezer.

BFUD-40-101 UPRIGHT FREEZER DOUBLE DOOR

Refrigerator: 225 L, Freezer: 225 L

LCD display

Audible/Visual alarm: High or low temperature alarm, system failure alarm

Upright type, ABS material interior, Outer body is painted steel board

60 mm extra thick insulation

The safety lock on the door, with door handle

4 Casters provided for easy handling

The international famous compressor and EBM fan motor

SPECIFICATIONS

Model	BFUD-40-101
Capacity	450 L
Temperature Range	Refrigerator: +2°C to +8°C, Freezer: -10°C~-40°C
Refrigerant	R507, CFC free
Internal Dimension	530x626x1356 mm
Overall Dimension	655x780x1845 mm
Drawers/Shelves	12
Weight (Net/Gross)	120 /140 kg
Average Power	450 W
Power Supply	220 V/110V, 50/60 Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com