



Product Image Coming Soon

## UPRIGHT FREEZER DOUBLE DOOR BFUD-40-502

## UPRIGHT FREEZER DOUBLE DOOR BFUD-40-502

High performance refrigerators featuring microprocessor control for temperature accuracy and ergonomic design for easy usage. Designed for space saving storage, it provides outstanding thermal performance along with safety features.

Used in Protection, Frozen Blood, Vaccines, Enzymes, Reagents, Research.

Also known as Double door refrigerated Freezer, Laboratory Combined Refrigerator and Freezer , Laboratory double door refrigerated freezer.

## BFUD-40-502 UPRIGHT FREEZER DOUBLE DOOR



Product Image Coming Soon

Vertical side opening door, free-standing door handle on the left side, easy to open  
4 casters for easy mobility; two of the casters have a brake function to prevent slippage

Dual LED digital display screen, can display the temperature data of the upper and lower rooms independently

Drawers are designed with PS materials, which are more durable and safer to store

## SPECIFICATIONS

Model	BFUD-40-502
Capacity	R:198L F:102L
Temperature Range	+2~8° & -10~-25°C
Alarm	High&low temp, power failure, control error, door open
External Dimension	650x620x1980 mm
Power	187~242V, 50/60Hz
Optional Parts	110V transformer; Datalogger



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)