



-25°C FREEZER CHEST BFUR-25-506

-25°C FREEZER CHEST BFUR-25-506

A perfect partner for your multi cold storage purposes in a contemporary lab set up. Microprocessor based temperature controls provide outstanding temperature uniformity with high quality compressor for energy efficiency. Safety features and durable interiors provide you with the most dependable performance.

Used in Laboratories, Routine Cold Storage, Biologicals, Pharmaceuticals, Clinical, Research.

Also known as Laboratory ULT freezer, ULT Freezer, Ultra-low Freezer.

BFUR-25-506 -25°C FREEZER CHEST



Undercounter Mini type Drawer x 4

SPECIFICATIONS

Model	BFUR-25-506
Capacity	100 L
Temperature Range	-10°C to -25°C
Alarm	High&low temp, power failure, control error, door open
External Dimension	540X540X880 mm
Power	187~242V, 50/60Hz
Optional Parts	110V transformer; Datalogger



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com