



## -40°C UPRIGHT FREEZER BFUR-40-308

## -40°C UPRIGHT FREEZER BFUR-40-308

Equipped with precise temperature controls, rugged interiors and safety features, our freezers provide you with the most dependable performance with low maintenance. A perfect partner for your multi cold storage purposes in a contemporary lab set up. Microprocessor based temperature controls provide outstanding temperature uniformity with high quality compressor for energy efficiency.

Used in Routine Cold Storage, Laboratories, Biologicals, Pharmaceuticals, Clinical, Research.

Also known as Laboratory ULT freezer, ULT Freezer, Ultra-low Freezer.

## BFUR-40-308 -40°C UPRIGHT FREEZER



Double Door, Drawer x 12, USB port

## SPECIFICATIONS

Model	BFUR-40-308
Capacity	358 L
Temperature Range	-20 to -40°C
Alarm	High&low temp, power failure, control error, door open, high ambient temperature
External Dimension	720×830×1870 mm
Power	187~242V, 50/60Hz
Optional Parts	110V transformer; Datalogger



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)