



## -40°C UPRIGHT FREEZER BFUR-40-309

## -40°C UPRIGHT FREEZER BFUR-40-309

Equipped with precise temperature controls, rugged interiors and safety features, our freezers provide you with the most dependable performance with low maintenance. A perfect partner for your multi cold storage purposes in a contemporary lab set up. Microprocessor based temperature controls provide outstanding temperature uniformity with high quality compressor for energy efficiency.

Used in Routine Cold Storage, Laboratories, Biologicals, Pharmaceuticals, Clinical, Research.  
Also known as Laboratory ULT freezer, ULT Freezer, Ultra-low Freezer.

## BFUR-40-309 -40°C UPRIGHT FREEZER



## SPECIFICATIONS

Model	BFUR-40-309
Capacity	305 L
Temperature Range	-20 to -40°C
Alarm	High&low temp, control error, power failure
External Dimension	1218x840x940 mm
Power	187~242V, 50/60Hz
Optional Parts	110V transformer; Datalogger



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)