



## ELECTRIC GYNECOLOGICAL TABLE BHBD-211

## ELECTRIC GYNECOLOGICAL TABLE BHBD-211

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.

Used in Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation.

Also known as Gyneacological table, Gyneacological bed, Examination table, Examination bed, Delivery bed, Gynaecological Delivery Bed.

## BHBD-211 ELECTRIC GYNECOLOGICAL TABLE



Whole stainless steel hidden auxiliary table with flexible extension and retraction

Equipped with hidden sewage basin

With foot control, can adjust the movements easily

Imported LINAK gear motor, without noise, have a long service life

## SPECIFICATIONS

Model	BHBD-211
Table Length	1450+50 mm
Table width	520+30 mm
Height Adjustment	(400 900)+50 mm
Trendelenburg/rev-trend	≥15°/≥25°
Standard Accessory	A pair of arm rods, a pair of leg support, a stainless steel sewage basin
Power Supply	220V 50Hz



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)