



## 3 CRANK MANUAL BED BHBD-510

## 3 CRANK MANUAL BED BHBD-510

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance

Used in Recuperating, Nursing.

Also known as Medical bed, hospital bed, Hospital cot, Electric ICU bed.

## BHBD-510 3 CRANK MANUAL BED



Steel framework, top class ABS & PP plastic

PP head & foot board and bed board, detachable and easy to clean 6 columns collapsible Aluminium alloy protective railing

4 pcs universal castor, dia. 125 mm, Individual brake

## SPECIFICATIONS

Model	BHBD-510
Function	3
Back-rest	0° - 65°
Knee Gatch	0° - 45°
HI - LOW	750 mm - 575 mm
Castor	360° swivel castors (Dia. 150 mm), with individual brake.
Material	Powder coated steel bed frame, top class ABS & PP plastic
Bed board dimension	1930x830 mm
Standard accessory	IV pole
Overall dimension	2160Lx970Wx575-750 mm
Patient surface	1930x830 mm
Height adjustment	750-575 mm
IV pole	Yes
Dimension	2160x960x530- 730 mm



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)