



## HOSPITAL TRANSFER STRETCHER TROLLEY BHBD-812

## HOSPITAL TRANSFER STRETCHER TROLLEY BHBD-812

Labstac Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

Used in Patient transfer, Patient trolley.

Also known as Medical stretcher, Ambulance stretcher, Emergency Stretcher, Patient transfer stretcher.

## BHBD-812 HOSPITAL TRANSFER STRETCHER TROLLEY



It is a type of manual transfer stretcher. Turn by deasil, the stretcher will be up, and turn by widdershins, the stretcher will be down;

Aluminium-Alloy protective rails, can be folded and hidden under the surface on no use;

Central-lock castors, diameter 150mm. Matched with direction castor, help the operator controlling the direction easily;

Rotate the crand and handle on bedboard, patient on the mattress can be moved parallelly to another bed or flat table

## SPECIFICATIONS

Model	BHBD-812
Height adjustment	920-610 mm
Net Weight	89 kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)