



## COVER SLIPER BHTP-503

## COVER SLIPER BHTP-503

Slide printer eliminates the need for writing and labeling by directly printing on the slides with precision. Used in Pathology, Cytology, Histology, Clinical Study Trials, Tracking of specimens. Also known as Slide printer, Cassettes printer.

## BHTP-503 COVER SLIPER



Covers of various brands and specifications. Color touch screen control interface, real-time dynamic prompts of operation instructions, simple and efficient operation, prevent misoperation. By adjusting the flow and speed of the sealing rubber, it can be compatible with imported and domestically produced multi-standard sealing glue, stepless adjustment and real-time response. Whether the dispensing needle is in real-time detection and alarm function at the working position, and a cleaning screen is designed at the storage position to effectively prevent the nozzle from drying (manual operation is required). Automatically monitor the quality of the coverslip during the working process. The fragments and bad pieces are automatically discarded into the waste cassette. The sensor has self-checking function to effectively improve work efficiency. Automatically carry out two rounds of slide scans during the work process, effectively reducing the leakage rate of the slides and improving the quality of work. The slide out box has the function of putting it at any time and taking it at any time, and the machine automatically performs the operation; the machine without the exit box automatically detects and has the graphic sound alarm function.

## SPECIFICATIONS

|                              |                           |
|------------------------------|---------------------------|
| Model                        | BHTP-503                  |
| Specimen slide throughput    | 200 ~ 400 pieces/hour     |
| Specimen slide holders       | 20 ~ 30 pieces/holder     |
| Output magazines             | 30 piece/magazine         |
| Cover slip magazine capacity | 120 ~ 160 pieces/magazine |
| Specimen slide size          | 26x75 mm                  |
| Cover slip size              | 24x50 mm                  |
| Weight                       | 70 kg                     |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)