



**Product Image Coming Soon**

## MINI COVER SLIPER BHTP-504

## MINI COVER SLIPER BHTP-504

Slide printer eliminates the need for writing and labeling by directly printing on the slides with precision. Used in Pathology, Cytology, Histology, Clinical Study Trials, Tracking of specimens. Also known as Slide printer, Cassettes printer.

## BHTP-504 MINI COVER SLIPER



Product Image Coming Soon

Covering machine is a fast curing device for sealed glass slides. This machine adopts advanced microcomputer technology, which makes the equipment more accurate, stable and anti-interference ability.

## SPECIFICATIONS

Model	BHTP-504
Slide seat	15
Speed	30 tablets/min
Rated voltage	AV 220 V $\pm$ 10 %
Power frequency	50 $\pm$ 1 Hz
Rated power	$\leq$ 100W



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)