



GENERAL PURPOSE INCUBATOR BIGP-201

GENERAL PURPOSE INCUBATOR BIGP-201

Biolab Scientific offers a range of high quality, efficient and reliable general purpose incubators that provide superior degree of temperature uniformity. Equipped with microprocessor PID controller, this series of incubators provide exceptional temperature control with added security of over temperature protection. These incubators serve the purpose for most biological analysis, research and general laboratory applications.

Used in Quality Control, Medicine Processing Industries, Research, Pharmaceuticals, Microbial Culturing.
Also known as Lab Incubator, Laboratory General Purpose Incubator.

BIGP-201 GENERAL PURPOSE INCUBATOR



Microprocessor controller (with Timer function)

With inner glass door for easy observation

Polished stainless-steel chamber

Independent temperature-limiting alarm system ensures experiments run safely.(Option)

Printer connector and RS485 connector are options which can connect printer and computer to record the parameters and the variations of temperature.(Option)

SPECIFICATIONS

Model	BIGP-201
Capacity	16 L
Temperature Range	RT +5-65°C
Ambient Temperature	+5-35°C
Temperature Uniformity	±1.5°C (@37°C)
Display Resolution	0.1°C /±0.5°C
Timer Range	0-5999 min
Internal Dimension	250Wx260Dx250H mm
Exterior Dimension	530Wx480Dx420H mm
Shelves	2
Power Consumption	200 W
Electrical Requirement	220V, 50Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com