



## GENERAL PURPOSE INCUBATOR BIGP-6903

# GENERAL PURPOSE INCUBATOR BIGP-6903

## BIGP-6903 GENERAL PURPOSE INCUBATOR



304 stainless steel, mirror polishing processing, easy to clean and maintain  
PID controller with over temperature alarm and timing function ensures precise and reliable control  
Air circulation system with specific air flow channel ensures a good temperature uniformity performance  
Silicon door gasket with long lifetime, and easy to change  
A damper adjustment in the front ensures the gas convection enough in working chamber

## SPECIFICATIONS

Model	BIGP-6903
Capacity	50 L
Electrical requirement	220 V 50 Hz
Temperature Range	RT +10 - 300°C
Display Resolution	0.1°C
Temperature Stability	±1°C
Temperature Uniformity	±3% (at 100°C)
Ambient Temperature	+5 - 40°C
Power Consumption	1100 W
Internal Dimension	420Wx395Dx350H mm
External Dimension	720Wx530Dx520H mm
Shelves	2 pcs
Timing Range	0~9999 min



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)