



GENERAL PURPOSE INCUBATOR BIGP-6904

GENERAL PURPOSE INCUBATOR BIGP-6904

BIGP-6904 GENERAL PURPOSE INCUBATOR



304 stainless steel, mirror polishing processing, easy to clean and maintain
PID controller with over temperature alarm and timing function ensures precise and reliable control
Air circulation system with specific air flow channel ensures a good temperature uniformity performance
Silicon door gasket with long lifetime, and easy to change
A damper adjustment in the front ensures the gas convection enough in working chamber

SPECIFICATIONS

Model	BIGP-6904
Capacity	80 L
Electrical requirement	220 V 50 Hz
Temperature Range	RT +10 - 300°C
Display Resolution	0.1°C
Temperature Stability	±1°C
Temperature Uniformity	±3% (at 100°C)
Ambient Temperature	+5 - 40°C
Power Consumption	1550 W
Internal Dimension	450Wx400Dx450H mm
External Dimension	740Wx530Dx630H mm
Shelves	2 pcs
Timing Range	0~9999 min



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com