



## GENERAL PURPOSE INCUBATOR BIGP-6906

# GENERAL PURPOSE INCUBATOR BIGP-6906

## BIGP-6906 GENERAL PURPOSE INCUBATOR



304 stainless steel, mirror polishing processing, easy to clean and maintain  
PID controller with over temperature alarm and timing function ensures precise and reliable control  
Air circulation system with specific air flow channel ensures a good temperature uniformity performance  
Silicon door gasket with long lifetime, and easy to change  
A damper adjustment in the front ensures the gas convection enough in working chamber

## SPECIFICATIONS

|                        |                   |
|------------------------|-------------------|
| Model                  | BIGP-6906         |
| Capacity               | 220 L             |
| Electrical requirement | 220 V 50 Hz       |
| Temperature Range      | RT +10 - 300°C    |
| Display Resolution     | 0.1°C             |
| Temperature Stability  | ±1°C              |
| Temperature Uniformity | ±3% (at 100°C)    |
| Ambient Temperature    | +5 - 40°C         |
| Power Consumption      | 2450 W            |
| Internal Dimension     | 600Wx500Dx750H mm |
| External Dimension     | 880Wx630Dx930H mm |
| Shelves                | 2 pcs             |
| Timing Range           | 0~9999 min        |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)