



GENERAL PURPOSE INCUBATOR BIGP-6906

GENERAL PURPOSE INCUBATOR BIGP-6906

BIGP-6906 GENERAL PURPOSE INCUBATOR



304 stainless steel, mirror polishing processing, easy to clean and maintain
PID controller with over temperature alarm and timing function ensures precise and reliable control
Air circulation system with specific air flow channel ensures a good temperature uniformity performance
Silicon door gasket with long lifetime, and easy to change
A damper adjustment in the front ensures the gas convection enough in working chamber

SPECIFICATIONS

Model	BIGP-6906
Capacity	220 L
Electrical requirement	220 V 50 Hz
Temperature Range	RT +10 - 300°C
Display Resolution	0.1°C
Temperature Stability	±1°C
Temperature Uniformity	±3% (at 100°C)
Ambient Temperature	+5 - 40°C
Power Consumption	2450 W
Internal Dimension	600Wx500Dx750H mm
External Dimension	880Wx630Dx930H mm
Shelves	2 pcs
Timing Range	0~9999 min



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com