



## GENERAL PURPOSE INCUBATOR BIGP-7602

# GENERAL PURPOSE INCUBATOR BIGP-7602

## BIGP-7602 GENERAL PURPOSE INCUBATOR



Microprocessor controller (with timing function)

With inner glass door for easy observation

Polished stainless-steel chamber

Independent temperature-limiting alarm system ensures experiments run safely

Printer connector and RS485 connector are options which can connect printer and computer to record the parameters and the variations of temperature

## SPECIFICATIONS

Model	BIGP-7602
Capacity	35 L
Electrical requirement	220 V 50 Hz
Temperature Range	RT+5 - 65°C
Display Resolution	0.1°C/±0.5°C
Ambient Temperature	+5 - 3°C
Temperature Uniformity	±1.5°C (at 37°C)
Power Consumption	200 W
Interior Dimension	340Wx320Dx320H mm
Exterior Dimension	620Wx490Dx490H mm
Shelves	2 pcs
Timing Range	0~5999 min
Remark	With observation window and without inner glass door



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)