



GENERAL PURPOSE INCUBATOR BIGP-7604

GENERAL PURPOSE INCUBATOR BIGP-7604

BIGP-7604 GENERAL PURPOSE INCUBATOR



Microprocessor controller (with timing function)

With inner glass door for easy observation

Polished stainless-steel chamber

Independent temperature-limiting alarm system ensures experiments run safely

Printer connector and RS485 connector are options which can connect printer and computer to record the parameters and the variations of temperature

SPECIFICATIONS

Model	BIGP-7604
Capacity	80 L
Electrical requirement	220 V 50 Hz
Temperature Range	RT+5 - 65°C
Display Resolution	0.1°C/±0.5°C
Ambient Temperature	+5 - 3°C
Temperature Uniformity	±1.5°C (at 37°C)
Power Consumption	400 W
Interior Dimension	500Wx400Dx400H mm
Exterior Dimension	780Wx530Dx560H mm
Shelves	2 pcs
Timing Range	0~5999 min
Remark	With observation window and without inner glass door



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com