



GENERAL PURPOSE INCUBATOR BIGP-7604

GENERAL PURPOSE INCUBATOR BIGP-7604

BIGP-7604 GENERAL PURPOSE INCUBATOR



Microprocessor controller (with timing function)

With inner glass door for easy observation

Polished stainless-steel chamber

Independent temperature-limiting alarm system ensures experiments run safely

Printer connector and RS485 connector are options which can connect printer and computer to record the parameters and the variations of temperature

SPECIFICATIONS

| | |
|------------------------|--|
| Model | BIGP-7604 |
| Capacity | 80 L |
| Electrical requirement | 220 V 50 Hz |
| Temperature Range | RT+5 - 65°C |
| Display Resolution | 0.1°C/±0.5°C |
| Ambient Temperature | +5 - 3°C |
| Temperature Uniformity | ±1.5°C (at 37°C) |
| Power Consumption | 400 W |
| Interior Dimension | 500Wx400Dx400H mm |
| Exterior Dimension | 780Wx530Dx560H mm |
| Shelves | 2 pcs |
| Timing Range | 0~5999 min |
| Remark | With observation window and without inner glass door |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com