



GENERAL PURPOSE INCUBATOR BIGP-7606

GENERAL PURPOSE INCUBATOR BIGP-7606

BIGP-7606 GENERAL PURPOSE INCUBATOR



- Microprocessor controller (with timing function)
- With inner glass door for easy observation
- Polished stainless-steel chamber
- Independent temperature-limiting alarm system ensures experiments run safely
- Printer connector and RS485 connector are options which can connect printer and computer to record the parameters and the variations of temperature

SPECIFICATIONS

| | |
|------------------------|--|
| Model | BIGP-7606 |
| Capacity | 270 L |
| Electrical requirement | 220 V 50 Hz |
| Temperature Range | RT+5 - 65°C |
| Display Resolution | 0.1°C/±0.5°C |
| Ambient Temperature | +5 - 3°C |
| Temperature Uniformity | ±1.5°C (at 37°C) |
| Power Consumption | 750 W |
| Interior Dimension | 600Wx600Dx750H mm |
| Exterior Dimension | 890Wx740Dx910H mm |
| Shelves | 2 pcs |
| Timing Range | 0~5999 min |
| Remark | With observation window and without inner glass door |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com