



## LABORATORY REFRIGERATOR BLAR-201

## LABORATORY REFRIGERATOR BLAR-201

Laboratory refrigerators supplied with digital controller, glass door, alarm system to meet the demanding requirements of laboratory research. These refrigerator has temperature ranging from +2°C to +8°C and all models are combined with two exact sensors and auto defrost.

Used in Industry, Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

Also known as Laboratory Medical Refrigerator.

## BLAR-201 LABORATORY REFRIGERATOR

Safety door lock

Auto defrost function

High ambient temperature and humidity area

Wide voltage from 192 V to 242 V

No condensation glass door

Functions: High / Low Temperature, High / Low Temperature Recorder, Power Failure, Low Battery, Door Ajar, Remote Alarm

Optional Accessories: Foaming door, Temperature Recorder

## SPECIFICATIONS

Model	BLAR-201
Capacity	60 L
Temperature Range	2~8°C
Cabinet Type	Desktop
Cooling Performance	4°C
Temperature Resolution	0.1°C
Refrigerant	R134a, CFC free ( R600a optional )
Defrost Method	Auto Defrost
Cooling Type	Forced air cooling
Material	Galvanized steel powder coating
Shelves	2
Climatic class	ST
Internal Size	426Wx388Dx435H mm
External Size	510Wx531Dx655H mm
Power	60 W
Power Supply	220V,50/ 60HZ, 110V,50/60HZ



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)