



## DESKTOP LAMINAR AIRFLOW BLDK-101

## DESKTOP LAMINAR AIRFLOW BLDK-101

Compact desktop design required for operations in an ultra clean, dust free environment. Small size saves precious laboratory space. Vertical and horizontal laminar flow cabinets available suit your needs.

Used in Medical, Pharmaceutical, Research, Manufacturing areas, Electronic setors, Industrial sectors, Laboratory, Healthcare.

Also known as Benchtop Laminar Air Flow, Laminar Flow Cabinet, Clean Bench, Laboratory Laminar Flow, Laboratory Laminar Air Flow, Tissue Culture Hood, Laminar Air Flow Cabinet, Minimal-Turbulence Air Flow, Laminar Flow Cabinet, Laminar Flow Hood, Laminar Flow.

## BLDK-101 DESKTOP LAMINAR AIRFLOW

Air Cleanliness: Class 100 @ $\geq 0.5 \mu\text{m}$   
Average Air Flow Velocity: 0.3-0.45 m/s  
Fluorescent Lamp: 8 W  
Noise Level:  $\leq 60$  dB(A)

## SPECIFICATIONS

Model	BLDK-101
Type	Horizontal Laminar Flow
Air Cleanliness	Class 100 @ $\geq 0.5 \mu\text{m}$
Average Air Flow Velocity	0.3-0.45 m/s
Fluorescent Lamp	8 W
Internal Size	600Wx495Dx500H mm
Overall Dimension	650x535x860 mm
Noise Level	$\leq 60$ dB(A)
Power	400 W
Power Supply	220 V, 50 Hz

## ACCESSORIES

Accessory Code	Name	Unit	Description
2800606006	Horizontal Laminar clean bench	1	
2800606007	Power cord	1	AC220V 60Hz
2800606008	Fuse - 5A	5	
2800606009	Fluorescent lamp- 14 W	1	
2800606010	Ultraviolet lamp - 14 W	1	



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)