



HORIZONTAL LAMINAR AIRFLOW BLHZ-201

HORIZONTAL LAMINAR AIRFLOW BLHZ-201

Compact design required for operations in an ultra clean, dust free environment. Ideal for laboratory applications where product protection is required. Small size saves precious laboratory space. Contains ultra thin filter including static pressure box without separator. Larger space permits working with laboratory equipments within the workspace. Used in Bioscience, Food processing, Pharmaceuticals, Aerospace, Medical research Laboratories, Hospitals, Research. Also known as Minimal-Turbulence Air Flow, Laboratory Laminar Flow, Clean Bench, Laboratory Laminar Air Flow, Laminar Flow Cabinet, Tissue Culture Hood, Laminar Air Flow Cabinet, Laminar Flow Cabinet, Laminar Flow Hood, Laminar Flow.

BLHZ-201 HORIZONTAL LAMINAR AIRFLOW



ETL certification, Satisfy North American market.
With memory function in case of power-failure.
Audio and visual alarm (Abnormal airflow velocity).

SPECIFICATIONS

Model	BLHZ-201
AirFlow Velocity	80~99fpm
HEPA Filter	One,99.995% efficiency at 0.3um
Pre filter	Polyster fibre, Washable
Work Surface Height	750 mm
Main Body Material	1.0 mm Cold-rolled steel with anti-bacteria powder coating.
Work Table Material	304 Stainless steel
Front Window	5 mm toughened glass, anti-UV.
LED Lamp	14Wx1
Display	LCD display
Waterproof Socket	One socket (Double type), America standard, Max.Power: 500W
Standard Accessory	Manually height adjustable base stand,LED lamp, wind speed sensor, waterproof socket , Stainless Steel LV.Bar with 10 hooks
Optional Accessory	UV lamp
Noise	268dB(A)
External Size	1100x810x (1720-1980)mm
Internal Size	1000Wx500Dx600H mm
Package Size	1240Wx1090Dx1250H mm
Gross Weight	190kg
Consumption	950W
Power Supply	103.5-126.5V, 60Hz, 10A

ACCESSORIES

Accessory Code	Name
6000610006	Fluorescent lamp
6000610007	UV lamp
6000610008	Base stand
6000610009	Waterproof socket



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com