



AUTOMATIC SEALING AND CAPPING MACHINE BLIH-102

AUTOMATIC SEALING AND CAPPING MACHINE BLIH-102

Automatic Sealing and Capping Machine is specifically designed to speed up the process of packaging of your product but also ensure the sterility of the product and that the product does not get contaminated

Used in Pharmaceuticals, Mobile inspection laboratories, Disease prevention and control centers, Food Industry, Beverage Industry.

Also known as Screw on capping machine.

BLIH-102 AUTOMATIC SEALING AND CAPPING MACHINE



One-button start
Simple and fast operation
Safe and labor-saving
Capping/sealing speed $\leq 4s$

SPECIFICATIONS

Model	BLIH-102
Throughput	8-96, suitable for 8-strip PCR tubes and 96-well PCR plates
Start Method	Button
Speed	Capping\sealing speed $\leq 4s$
Function	8-strip PCR tubes capping, 96-well PCR plate capping, PCR plate sealing film
Adapter	2 types, 8-strip PCR tubes and 96-well PCR plates, can be customized
Pressure	Peak pressure: 300kg
External Size	277x194x230 mm
Package Size	300x260x340 mm
Weight	Net Weight- 9.2kg Gross Weight- 12kg
Power Supply	100-240V 50/60Hz



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com