



MULTI-PARAMETER ANALYZER BMET-1100

MULTI-PARAMETER ANALYZER BMET-1100

Multi-parameter Analyzer simultaneously measures Ion, Conductivity, Dissolved Oxygen, Temperature that otherwise requires four different meters. With the appropriate probes installed, it can also simultaneously measure several or one parameters.

BMET-1100 MULTI-PARAMETER ANALYZER



SPECIFICATIONS

Model	BMET-1100
DC Voltage	
600 mV/6 V/60 V/600 V/1000 V	± (0.5%+3)
AC Voltage	
6 V/60 V	± (0.8%+3)
600 V/750 V	± (1.0%+10)
DC Current	
60 µA/60 mA/600 mA	± (0.8%+3)
20 A	± (1.2%+3)
AC Current	
60 mA/600 mA	± (1.0%+3)
20 A	± (1.5%+3)
Resistance(Ω)	
600 /6 k/6 0k/600 k/6 M	± (0.8%+3)
60 MΩ	± (1.2%+30)
Capacitance	
6 nF	± (4.0 % +30)
60 nF/600 nF/6 µF/60 uF/600 µF/6 mF ±(4.0%+3)	± (4.0 % +3)
100 mF	± (5.0 % +3)
Frequency	
60/600/6 k/60 k/600 k/6 M/10 MHz	± (1.0%+3)
Duty	
5 % ~95 %	± (2.0 % +3)
hFE	
0 ~1000	± (1.0 % +3)



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com