



MULTI-PARAMETER ANALYZER BMET-1106

MULTI-PARAMETER ANALYZER BMET-1106

Multi-parameter Analyzer simultaneously measures Ion, Conductivity, Dissolved Oxygen, Temperature that otherwise requires four different meters. With the appropriate probes installed, it can also simultaneously measure several or one parameters.

BMET-1106 MULTI-PARAMETER ANALYZER



SPECIFICATIONS

Model	BMET-1106
DC Voltage	
200 mV/2 V/20 V/200 V	$\pm (0.5\%+2)$
600 V	$\pm (0.8\%+2)$
AC Voltage	
2 V/20 V/200 V	$\pm (1.2\%+10)$
DC Current	
200 μ A/2 mA/20 mA	$\pm (1.0\%+2)$
200 mA	$\pm (1.5\%+2)$
10 A	$\pm (3.0\%+2)$
Resistance	
200 Ω /2 k Ω /20 k Ω /200 k Ω	$\pm (0.8\%+2)$
2 M Ω	$\pm (1.0\%+2)$
Display	2000 Counts
Data Hold	Yes
Diode Test	Yes
Continuity	Yes
Power Supply	1x9V 6F22 Batteries
Weight	Approx. 169.5/173/174.5 g
Size	148x74x42 mm
Safety Standard/Rating	EN61010-1,EN61010-2-033, EN61326,CAT.III 600 V



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com