



MULTI-PARAMETER ANALYZER BMET-1107

MULTI-PARAMETER ANALYZER BMET-1107

Multi-parameter Analyzer simultaneously measures Ion, Conductivity, Dissolved Oxygen, Temperature that otherwise requires four different meters. With the appropriate probes installed, it can also simultaneously measure several or one parameters.

BMET-1107 MULTI-PARAMETER ANALYZER



SPECIFICATIONS

Model	BMET-1107
DC Voltage	
200 mV/2 V/20 V/200 V	$\pm (0.5\%+2)$
600 V	$\pm (0.8\%+2)$
AC Voltage	
2 V/20 V/200 V	$\pm (1.0\%+3)$
600 V	$\pm (1.2\%+3)$
DC Current	
200 μ A/2 mA/20 mA/200 mA	$\pm (1.5\%+3)$
AC Current	
200 μ A/2 mA/20 mA/200 mA	$\pm (1.5\%+4)$
DC	
10 A	$\pm (3.0\%+3)$
AC	
10A	$\pm (3.0\%+3)$
Resistance	
200 Ω /2 k Ω /20 k Ω /200 k Ω /2 M Ω	$\pm (0.8\%+2)$
20M Ω	$\pm (1.0 \% + 2)$
Temperature	
- 20 °C ~ 1000 °C	$\pm (3.0 \% + 3)$
0 °F ~ 1800 °F	$\pm (3.0 \% + 3)$
Display	2000 Counts
Data Hold	Yes
Continuity	Yes
Diode Test	Yes
Maximum Value Hold	Yes

Power Supply	1x9 V 6F22 Battery
Weight	Approx. 172.5/173 g
Size	148x74x42 mm
Safety Standard/Rating	ENG1010-1,ENG1010-2-033, EN61326,CAT III 600 V



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com