



MULTI-PARAMETER ANALYZER BMET-1112

MULTI-PARAMETER ANALYZER BMET-1112

Multi-parameter Analyzer simultaneously measures Ion, Conductivity, Dissolved Oxygen, Temperature that otherwise requires four different meters. With the appropriate probes installed, it can also simultaneously measure several or one parameters.

BMET-1112 MULTI-PARAMETER ANALYZER



SPECIFICATIONS

Model	BMET-1112
DC Voltage	
400 mV/4 V/40 V/400 V	$\pm (0.8\%+2)$
600 V	$\pm (0.1\%+2)$
AC Voltage	
40 V/400 V/600 V	$\pm (1.2\%+10)$
DC Current	
400 μ A/4 mA/40 mA/400 mA	$\pm (0.1\%+2)$
400 mA	$\pm (1.5\%+2)$
10 A	$\pm (3.0\%+2)$
Resistance	
400 Ω /4 k Ω /40 k Ω /400 k Ω	$\pm (0.8\%+3)$
Display	4000
Manual Rang	Yes
Bargraph Display	Yes
Data Hold	Yes
Diode Test	Yes
Continuity	Yes
NCV Non-Contact Voltage Detection	Yes
Back Light	Yes
Flashlight	Yes
Power Supply	2 x 1.5VAAA
Weight	268.5
Size	166x78x48 mm
Safety Standard/Rating	EN61010-1,EN61010-2-033,EN61326,CAT.III 600 V



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com