



MAGNETIC STIRRER BMGS-102

MAGNETIC STIRRER BMGS-102

Magnetic stirrer series offers many work modes, Strong magnetism, 30° bevel control panel, Brushless DC motor for convenient operations. Automatic reverse rotation for better-mixing result. Slip-proof and stable. Programmable memory can store nine common set S1 ~ S9, Meet the different kinds of experimental demands.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Magnetic Stirrer, Digital Magnetic Stirrer.

BMGS-102 MAGNETIC STIRRER

Single-channel

Specially designed for cell culture that needs low speed stirring

Programmable, make the operation mode automatically

Digital LED displays the actual speed

Many work modes, including the constant speed stirring, cycle stirring, two speed stirring and combination cycle and two speed stirring

SPECIFICATIONS

Model	BMGS-102
Capacity	3 L
Stirring Positions	1
Speed Range	5-200 rpm
Dimension	205x250x90 mm
Power	15 W
Power Supply	100-240V, 50/60Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Description	Capacity
3400807006	Spinner Flask	Standard; Flat Bottom; Complete	500 ml
3400807007	Spinner Flask	Standard; Flat Bottom; Complete	3000 ml



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com