



## BIOLOGICAL MICROSCOPE BMIC-302

## BIOLOGICAL MICROSCOPE BMIC-302

Digital Biological Microscope is an optical microscope that is primarily designed to observe cells, tissues as well as other biological specimens. Magnification range anywhere from 10x - 1,000x or more by attaching multiple objective lenses.

Used in University, Medical, Laboratory, Research, Pharmaceutical.

## BMIC-302 BIOLOGICAL MICROSCOPE



The model has a sliding trinocular viewing head at 45° which is 360° rotatable. Valid mould proof design. It can be used under the high temperature and high humidity environment.

## SPECIFICATIONS

Model	BMIC-302
Viewing Head	
Sliding binocular head at 45°, 360° rotatable	NA
Sliding trinocular head at 45°, 360° rotatable	Standard
Eyepiece	
Eyepiece WF10 x/ 18 with scale of crosshair	Standard
Eyepiece WF 10 X/18	Standard
Objective	
Strain free achromatic objective 4X , 10X , 40X	Standard
Nosepiece	
Toward quadruple nosepiece, center adjustable	Standard
Analyzer	
Rotatable analyzer with gradation 0°- 90°	Standard
Bertrand Lens	
Bertrand lens, sliding in/ out of optical path	Standard
Optical Compensator	
$\lambda$ slip(first class red)	Standard
$\frac{1}{4} \lambda$ slip	Standard
Quartz wedge	Standard
Revolving Round Stage	
Diameter $\Phi$ 160 mm, graduated in 1° increments, Minimum resolution 6' when using vernier scale	Standard
Condenser	
Abbe condenser with Iris diaphragm & filter	Standard
Focusing	
Coaxial coarse & fine adjustment, range 28 mm, fine division 0.002 mm	Standard

Polarizer	
Sliding in/out of optical path, located on the top of collector	Standard
Illumination	
6 V / 20 W halogen lamp, brightness adjustable	Standard
Power Supply	
AC 110-220 V , 50/60 Hz	Standard
Package Size	
327x273x425 mm	Standard
Gross Weight	
7.5kg	Standard



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)