



MULTISTATION STIRRER BMSS-104

MULTISTATION STIRRER BMSS-104

Powerful Multi-Station magnetic stirrers with 4 to 16 stirring positions to minimize bench space. Uniform agitation reaction conditions are achieved by synchronized operation. Speed ranges from 0-1100 rpm. Individually controlled stirring for consistency with various samples.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Multistation Magnetic Stirrer, Digital Magnetic Stirrer.

BMSS-104 MULTISTATION STIRRER



Speed range of 0-1100 rpm

High performance magnetic stirrer with 10 stirring positions benefit laboratory space

Individually controlled stirring for consistency with various samples

Brushless DC motor, maintenance free and provides constant speed

Stainless steel plate with silicone film, anti-slip and anti-corrosion

SPECIFICATIONS

Model	BMSS-104
Capacity	400 ml
Stirring Positions	10
Speed Range	0~1100 rpm
Temperature	Permissible Ambient Temperature-5-40°C
Permissible Relative Humidity	0.8
Work Plate Size	180x450 mm
Platform Material	Stainless steel with silicon film
Speed Display	Scale
Motor Type	Brushless DC motor
Dimension (WxDxH)	182x552x65 mm
Weight	7 kg
Power Supply	AC110/220V, 50/60Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com