



MULTISTATION STIRRER BMSS-301

MULTISTATION STIRRER BMSS-301

Powerful Multi-Station magnetic stirrers with 4 to 16 stirring positions to minimize bench space. Uniform agitation reaction conditions are achieved by synchronized operation. Speed ranges from 0-1100 rpm. Individually controlled stirring for consistency with various samples.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Multistation Magnetic Stirrer, Digital Magnetic Stirrer.

BMSS-301 MULTISTATION STIRRER



Ultra-light and ultra-thin, novel appearance, compact size, saving experiment space
The controller and the agitator are separated. The number and combination of stirrer can be adjusted arbitrarily, up to 6 stirrer

Adopt the connection method of USB interface, easy to plug and unplug

Each stirrer has its own speed and time, which makes it easier to compare experiments

Touch keys and sliding touch adjustment, durable, novel and fashionable

Advanced non-wearing magnetic coil technology, the operation is quieter and more stable, and the service life is longer

The LCD screen displays the parameters in real time, the LED indicator shows the working status, and the instrument running status is more intuitive

SPECIFICATIONS

Model	BMSS-301
Controller	
Max. stirring volume (H2O)	/
Motor output power	/
Speed setting (running) range	20-1500 rpm (increment 10)
Timing Range	0-9 hours(0 means infinity, increment 5min)
Display method	LCD
Working plate material	/
Working plate size	/
Dimensions	165x105x45 mm
Weight	0.48 Kg
Power adapter input	AC100-240V 50/60Hz 60VA
Power adapter output	DC 24V, 2A
Stirrer	
Max. stirring volume (H2O)	1 L
Motor output power	4 W
Speed setting (running) range	20-1500 rpm
Timing range	/
Display method	/
Working plate material	PVC
Working plate size	Φ120mm
Dimensions	120x120x19 mm
Weight	0.18 Kg
Power adapter input	/

Power adapter output	/
----------------------	---

OPTIONAL ACCESSORIES

Accessory Code	Name	Description	Waterproof
3401107006	Controller	Controller	No
3401107007	Stirrer Plate	20 ~1500 rpm	No
3401107008	Stirrer Plate	20 ~1000 rpm	Yes



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com