



GENERAL PURPOSE OVEN BOGP-102

GENERAL PURPOSE OVEN BOGP-102

Designed for continuous drying operations , general purpose ovens provide temperature homogeneity regulated by microcomputer control system. Stainless steel interiors with easy to clean surface while silicon rubber sealing prevents heat loss and gives stability during operation.

Used in Drying, Roasting, Laboratory, Research, Engineering, Industry, Agricultural, Biology.

Also known as Laboratory Oven, Laboratory general Purpose Oven.

BOGP-102 GENERAL PURPOSE OVEN



Micro-computer intelligent control system with LED display

Forced air convection maintains temperature uniformity

High quality stainless steel inner chamber with electrostatic spray treated exterior

Large observation window

Easy to clean semicircular corners with adjustable shelves

Silicon rubber sealing ring prevents heat loss

Over temperature alarm with timer

SPECIFICATIONS

Model	BOGP-102
Capacity	70 L
Temperature Range	RT+5-250°C
Temperature Accuracy	±1°C
Temperature Uniformity	±4.5°C at 250°C
Temperature Resolution	±1°C
Timer Range	0-999 min
Internal Dimension	420Wx450Dx350H mm
Exterior Dimension	630Wx740Dx550H mm
Package Size	750x820x630 mm
Shelves	2
Weight	44.5/47.5 kg
Power	1200 W
Power Supply	220V, 60Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com