



NATURAL CONVECTION OVEN BONC-101

NATURAL CONVECTION OVEN BONC-101

With easy to clean interiors and natural convection circulation system for efficient heating of chamber, this oven can be relied on for general laboratory applications.

Used in Drying, Roasting, Laboratory, Research, Engineering, Industry, Agricultural, Biology.

Also known as Laboratory Oven, Laboratory natural Convection Oven.

BONC-101 NATURAL CONVECTION OVEN



304 stainless steel, mirror polishing processing, easy to clean and maintain

Natural convection with low noise

PID controller with over temperature alarm and Timer function ensures precise and reliable control

SPECIFICATIONS

Model	BONC-101
Capacity	27 L
Temperature Range	RT+10-200°C
Ambient Temperature	+5-35°C
Temperature Stability	±1°C
Display Resolution	0.1°C
Internal Dimension	320Wx300Dx355H mm
Exterior Dimension	460Wx520Dx660H mm
Shelves	2
Power consumption	850 W
Electrical Requirement	AC 220V, 50Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com