



NATURAL CONVECTION OVEN BONC-308

NATURAL CONVECTION OVEN BONC-308

With easy to clean interiors and natural convection circulation system for efficient heating of chamber, this oven can be relied on for general laboratory applications.

Used in Biochemical, Pharmaceutical, Healthcare, Agricultural, Environmental protection, disinfection, Sterilization, Drying, Roasting, Laboratory, Research, engineering, Industry,.

Also known as Vacuum Oven, Laboratory Vacuum Oven.

BONC-308 NATURAL CONVECTION OVEN



Stainless steel inner chamber without scroll rolled

Silicone rubber door seal

Strong sealing

Nice design

SPECIFICATIONS

Model	BONC-308
Capacity	25 L
Temperature Accuracy	$\pm 1^{\circ}\text{C}$
Temperature Fluctuation	$\pm 2\%$
Vacuum Degree	$< 133 \text{ Pa}$
Timer Range	0-999 min
Internal Dimension	300Wx300Dx280H mm
Shelves/Trays	2
Power	900 W
Power Supply	220 V, 50 Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com