



## HOLE-TYPE OR LAMP BOPL-204

## HOLE-TYPE OR LAMP BOPL-204

Hole type operation light is suitable for small operating room. Labstac lights uses multiple different sources and by focusing knob, the beam is illuminated from different surgery locations which reduces the shadow by medical staff head, hands and other parts caused. The illumination adjusted according to actual needs.

Used in Medical unit, Surgery, Operation theatre.

Also known as Hole type surgical Light, Hole type OT light.

## BOPL-204 HOLE-TYPE OR LAMP



## SPECIFICATIONS

Model	BOPL-204
Illumination	$\geq 160,000\text{LUX}$
Color temperature	$4500\text{K} \pm 500\text{K}$
Color rendering Index	$100 \geq \text{Ra} \geq 85$
Depth of light beam	$> 500\text{ cm}$
The Value of total irradiance	$< 6\text{ mW/m}^2\text{lx}$
Total irradiance	$< 600\text{ W/m}^2$
Temperature rise on the head of the patient under surgical operation	$\leq 2^\circ\text{C}$
Rated Voltage of Bulb	24V
Rated Power of bulb	300 W
Input power	350 VA
Power Supply	$220\text{V} \pm 10\% 50\text{Hz}$



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)