



## LED OR SHADOWLESS LAMP BOPL-603

## LED OR SHADOWLESS LAMP BOPL-603

Challenging tasks require special solutions. The innovative Labstac LED operation light give you benefits of cutting-edge LED technology in a uniquely effective manner. We focus on providing LED lights that are more color-accurate, more flexible, more economical. We offer innovative light and system solutions that can be retrofitted easily and without special preparation.

Used in Dental implant surgery, Medical unit, Surgery, Operation theatre.

Also known as Surgical LED Light, LED operation theatre light, LED OT light.

## BOPL-603 LED OR SHADOWLESS LAMP



Adopt imported LED cold light source, white the ultra-thin optical lens

Adopt digital means for stepless regulation of LED brightness

Adjust the diameter of light spot and color temperature

With Endoscopy model

Removable handle sheath, easy to clean and disinfect

Nano coating of heat sink has created an excellent heat dissipation effect.

Imported balance arm, which is flexible up and down

## SPECIFICATIONS

Model	BOPL-603
Illumination	40000~160000 LUX (at 1m distance)
Color temperature	6700>TC>3000K Endo model- 10000 lux
Color rendering Index	100≥Ra≥85
Working distance	70-140 cm
Diameter of spot D50	110+35 mm
Light field diameter	10-28 cm (at 1m distance)
Heat to light ratio	<6mW/m <sup>2</sup> lx
D10	180±40 mm
Number of LED bulbs	45
Average bulb life (hours)	50000 hours
Max power consumption	90 VA



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)