



## MANUAL OR TABLE BOPT-201

## MANUAL OR TABLE BOPT-201

Manual Operation table are suitable for various operation tables conducted by medical unit. The OR table and accessories are ergonomic and can be used across multiple disciplines. Various actions are controlled by hydraulic pneumatic system. 360° rotation helps every surgical requirement.

Used in Medical unit, Surgery, Operation theatre, cardiovascular, pediatric, gynecology, abdominal, gall bladder, plastic surgery and more.

Also known as Manual surgical table, Manual OR Table, Manual operating table.

## BOPT-201 MANUAL OR TABLE



The super low level design can enable doctor to conduct operation sitting by operating table

Adopt electric hydraulic double cylinder elevating with stable performance

Import rotor and Electromagnetic valve

Imported sealing ring and is durable

The base seat and upright column covering are made of high quality 304 stainless steel which is anti-corrosive and can be easily cleaned

Electric braking control which is easily for moving operating table

## SPECIFICATIONS

Model	BOPT-201
Table Length	2100±50 mm
Table Width	480±20 mm
Height Adjustment	(800~1050) ±20 mm
Head board upwards	≥45°
Head board downwards	≥90°
Backboard upwards	≥80°
Backboard downwards	≥30°
Waist board lifting	100+20 mm
Leg board downwards	≥90°
Back Board Upwards	≥80°
Back Board Downwards	≥30°
Standard Accessories	A Pair of Shoulder/Waist Support, A pair of Arm Support, A pair of leg Plate, A screen



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)