



## ORTHOPEDICS TRACTOR RACK BOPT-401

## ORTHOPEDICS TRACTOR RACK BOPT-401

Biolab's Orthopedics tractor rack is mainly applied to orthopedics surgical by using this tractor rack the patient's leg can be performed side traction, prone traction and open-up position of the body during the operation.

Used in Surgery, Leg operation, Medical unit.

Also known as Operating Table.

## BOPT-401 ORTHOPEDICS TRACTOR RACK



Adopt carbon fiber material, the prespective imaging is clear, light and strong

The plug in design is convient and quick to connect with table

Equipped with a traction frame accessory car, which is convient to move and save space

## SPECIFICATIONS

|                                   |             |
|-----------------------------------|-------------|
| Model                             | BOPT-401    |
| Table Length                      | 1420±50 mm  |
| Table Width                       | 765±50 mm   |
| Traction device horizontal rotary | 0-180°      |
| Traction distance                 | (0-200)+10% |
| Traction Device height adjustment | >400 mm     |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)