



## VACUUM OVEN BOVA-7206

# VACUUM OVEN BOVA-7206

## BOVA-7206 VACUUM OVEN

304 stainless steel, mirror polishing processing, easy to clean and maintain

PID controller with over temperature alarm and timing function ensures precise and reliable control, also save more than 40% heating time

A big dual layer tempered glass on the door provides good observation

Door adjustment system with silicon door gasket ensures better vacuity

## SPECIFICATIONS

Model	BOVA-7206
Capacity	50 L
Electrical requirement	AC 220 V 50 HZ
Controller	LED display
Power Consumption	1450 W
Temperature range	RT-10 - 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	415Wx370Dx345H mm
Exterior Dimension	730Wx560Dx550H mm
Shelves	3 pcs
Chamber Material	304 Stainless steel



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)