



VACUUM OVEN BOVA-7219

VACUUM OVEN BOVA-7219

BOVA-7219 VACUUM OVEN

304 stainless steel, mirror polishing processing, easy to clean and maintain

PID controller with over temperature alarm and timing function ensures precise and reliable control, also save more than 40% heating time

A big dual layer tempered glass on the door provides good observation

Door adjustment system with silicon door gasket ensures better vacuity

SPECIFICATIONS

Model	BOVA-7219
Capacity	30 L
Electrical requirement	AC 220 V 50 HZ
Controller	LED display
Power Consumption	300 W
Temperature range	RT-10 - 65°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	320Wx320Dx300H mm
Exterior Dimension	610Wx510Dx490H mm
Shelves	1 pcs
Chamber Material	304 Stainless steel



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com