



## VACUUM OVEN BOVA-7304

# VACUUM OVEN BOVA-7304

## BOVA-7304 VACUUM OVEN



Microprocessor controller

Dual layer tempered glass door for clear observation

Minimum heating time 50% less than traditional vacuum oven

Interior chamber made from stainless steel with a mirror polished finish making for easy cleaning

Standard with vacuum pump

## SPECIFICATIONS

Model	BOVA-7304
Capacity	215 L
Electrical requirement	AC 220 V 50 HZ
Controller	LED display
Power Consumption	2100 W
Temperature range	RT+10 - 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	560Wx600Dx640H mm
Exterior Dimension	720Wx8200Dx1750H mm
Shelves	6 pcs
Chamber Material	304 Stainless steel



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)