



## VACUUM OVEN BOVA-7306

# VACUUM OVEN BOVA-7306

## BOVA-7306 VACUUM OVEN



Microprocessor controller

Dual layer tempered glass door for clear observation

Minimum heating time 50% less than traditional vacuum oven

Interior chamber made from stainless steel with a mirror polished finish making for easy cleaning

Standard with vacuum pump

## SPECIFICATIONS

Model	BOVA-7306
Capacity	90 L
Electrical requirement	AC 220 V 50 HZ
Controller	LED display
Power Consumption	1350 W
Temperature range	RT+10 - 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	450Wx450Dx450H mm
Exterior Dimension	610Wx590Dx1350H mm
Shelves	4 pcs
Chamber Material	304 Stainless steel



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)