



PIPETTE PUMP BPIC-201

PIPETTE PUMP BPIC-201

We are engaged in offering a wide range of Pipette Pumps that provide safe and accurate pipetting and dispensing of liquids. It can be rotated to dispense and draw up the liquid. Available in variety of sizes and shapes.

Used in Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research.

Also known as Laboratory Pipette Pump.

BPIC-201 PIPETTE PUMP



Economical
Highly gas tightness
Environment- friendly materials
Lightweight and dexterous for operation
Durable
Recalibration Mechanism

SPECIFICATIONS

Model	BPIC-201
Volume	2ml
Color	Bule + White
Unit	1000



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com