



## PIPETTE PUMP BPIC-202

## PIPETTE PUMP BPIC-202

We are engaged in offering a wide range of Pipette Pumps that provide safe and accurate pipetting and dispensing of liquids. It can be rotated to dispense and draw up the liquid. Available in variety of sizes and shapes.

Used in Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research.

Also known as Laboratory Pipette Pump.

## BPIC-202 PIPETTE PUMP



Economical  
Highly gas tightness  
Environment- friendly materials  
Lightweight and dexterous for operation  
Durable  
Recalibration Mechanism

## SPECIFICATIONS

Model	BPIC-202
Volume	10 ml
Color	Green + White
Unit	1000



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)