

SEROLOGICAL PIPETTE BPIC-607


SEROLOGICAL PIPETTE BPIC-607

Manufactured for accurate and rapid fluid dispensing. It is sterile, pyrogen-free, cotton plugged, and can be used for a variety of liquid handling applications. The reservoir eliminates the need for multiple reagent containers and saves bench top space.

Used in Plasmid Purification, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

Also known as Laboratory Serological Pipette.

BPIC-607 SEROLOGICAL PIPETTE



Easy to use

Medical grade original Polystyrene

RNase- free, DNase- free and Non-pyrogenic

Seven types of volume 1ml, 2ml, 5ml, 10ml, 25ml, 50ml and 100ml

Bidirectional graduation improve accuracy

Seven types of color code on tips make pipette

Sterilized by Gamma irradiation

SPECIFICATIONS

Model	BPIC-607
Volume	10 ml
Precision	2%
Negative Graduation (ml)	-3
Graduation	10-Jan
Color	Red
Sterile	Gamma
Package	Individual
Pcs/bag	25 pcs/bag
Pcs/ctn	4000 pcs/carton
Unit	1



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com