

# SEROLOGICAL PIPETTE BPIC-610


## SEROLOGICAL PIPETTE BPIC-610

Manufactured for accurate and rapid fluid dispensing. It is sterile, pyrogen-free, cotton plugged, and can be used for a variety of liquid handling applications. The reservoir eliminates the need for multiple reagent containers and saves bench top space.

Used in Plasmid Purification, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

Also known as Laboratory Serological Pipette.

## BPIC-610 SEROLOGICAL PIPETTE



Easy to use

Medical grade original Polystyrene

RNase- free, DNase- free and Non-pyrogenic

Seven types of volume 1ml, 2ml, 5ml, 10ml, 25ml, 50ml and 100ml

Bidirectional graduation improve accuracy

Seven types of color code on tips make pipette

Sterilized by Gamma irradiation

## SPECIFICATIONS

Model	BPIC-610
Volume	25 ml
Precision	2%
Negative Graduation (ml)	-3
Graduation	05-Jan
Color	Purple
Sterile	Gamma
Package	Bulk pack
Pcs/bag	25 pcs/bag
Pcs/ctn	2100 pcs/carton
Unit	1



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)