



TIPS BOX FOR 10µL TIPS BPIC-701

TIPS BOX FOR 10?L TIPS BPIC-701

It offers quick and easy access to pipette tips, saving valuable time. Space-saving and environmentally compatible. Commonly used in liquid handling and molecular biology.

Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

BPIC-701 TIPS BOX FOR 10?L TIPS



Made from medical grade Polypropylene

Slippy surface with excellent design

DNase-free

Rnase-free

Non-pyrogenic

Filter and non-filter for choice

Bag and rack are available to pack

SPECIFICATIONS

| Model | BPIC-701 |
|--------|----------|
| Volume | 10 µl |
| Unit | 20 |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com